


Group Varian

Insulin Pump UI Proposal

Patric Weber
Frank Kirchhoff
Hari Bajekal
Stephan Renggli



Dashboard

The Dashboard shows a quick overview of the most important data and acts as the main menu of our app.


Basal dose in U/h

Bolus dose in Units

Blood glucose Level
in mmol/L


Profile selector

Status & Settings



Basal

The Basal screen displays the currently set basal value and a histogram of past values.



Bolus

The Bolus screen allows you to input the carb count of a meal and administer the correct amount of Insulin.

Carb input in kcal

Administer button

Number input field

Dashboard Bolus

Enter carb count of your meal.

800 kcal

Administer Insulin

1 2 3


4 5 6

7 8 9

0


Blood Glucose

The Blood Glucose screendisplays the curren blood glucose level and a histogram of the last 24h.



Current blood glucose level


Blood glucose level histogram displaying the last 24h



Profiles


The Profiles screen allows you to set up and select multiple profiles for different situations in your life.

Multiple editable profiles
that control basal, bolus
and other values




Status

The Status screen contains information about your device such as battery level, connection status and the condition of the Insulin.



Battery,
Bluetooth connection and
Insulin condition



Settings

The Settings screen will help you set up your personal values such as a Medical ID, medication and other preferences.

Medical ID

Various settings
and preferences

